

PY-130 DeBlistering Machine



**EQUIPMENT
TECHNICAL PARAMETER**

南京翰洋智能设备有限公司
NANJING HANYOO PHARMATECH CO., LTD.



PY-130 DeBlistering Machine



INTRODUCTION

The deblistering machine is an automatic tablet or capsule deblistering machine for recovering drugs in the defective aluminum plastic blister packaging drug board. It is in the 304 grade stainless steel housing, micro motor, a reducer, roller combination, guide combination block rubber plate sequential connection, and equipped with electrical brake switch and break lock

switch. The micro-motor and a first-stage reducer drive a set of rollers that simulate the manual peeling action of the tablets, peel the blister aluminum foil tablets of the whole board from the defective packaging board, and collect the pills and packaging separately through the guide plate and the grain rubber plate.

FEATURES

1. Good appearance, simple and generous, not complicated, easy to clean, humanized design, easy to use and safe.
2. 304 Stainless steel, excellent material; Using pure PU material roller, wear-resistant, will not drop scraps, maximize the protection of drugs from any secondary pollution
3. The equipment has high precision, good use effect, capsules and tablets are completely out of the aluminum-plastic plate, will not crush the capsules and tablets, and the service life of the equipment is long.
4. Reasonable design, attention to details, the whole design process is from the actual user's point of view, safe operation, easy to move, different sizes, arrangement of the plate adjustment is convenient, no need to replace the mold. Maximum close to the actual use of workshop needs.
5. The products have been recognized by European and American customers, and are highly praised by domestic users, with good effect, small and flexible, convenient and practical. European and American quality, mid-range price, cost-effective.

TECHNICAL PARAMETER

MODEL	PY-130
Working Mode	Full-automatic
Counting Mode	Automatic
Work Efficiency	100-120 boards/min
Scope Of Application	Capsules, tablets, soft capsules, candy and other all aluminum packaging, aluminum packaging, aluminum packaging
Application Plate	The grain should be arranged horizontally and vertically
Change Mold	No mold required
Driving Part	Linear Guide rail - Taiwan HIWIN
Motor And Governor	Taizhou TAILI Taili
Power Requirements	220V AC:50-60HZ
Rated Power	180W (Adjustable Speed)
Equipment Weight	38KG
Overall Dimension	400*400*1100mm
Equipment Material	Stainless steel304

The banner features a dark blue background with several white line-art icons related to medicine and healthcare. These include a blister pack of pills, a medicine bottle with a cross, a pill being split, a pill bottle, two pills, a syringe, and another blister pack.

HANYOO

翰洋科技

Add:No.380 Sanqiaogeng,Gaochun District,Nanjing,China

Tel:+86-025-5621 6295

[Http://www.hanyoopm.com](http://www.hanyoopm.com)